*Guida didattica / Didactic guide*  
*Guía didáctica*

Ver. 1.0.1

---

**Cod. 2056**

---

**OPTIKA S.R.L.**

VIA RIGLA, 30 – 24010 PONTERANICA (BERGAMO) – ITALY  
Tel. +39 035 571392 - Fax +39 035 571435

[www.optikascience.com](http://www.optikascience.com)

[info@optikascience.com](mailto:info@optikascience.com)

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## WATER CALORIMETER

This calorimeter has a capacity of approximately 1000 ml. It is supplied with a thermometer and a stirrer.

### 1) DETERMINING THE WATER EQUIVALENT OF THE CALORIMETER

During heat exchanges in the calorimeter you have to consider that part of the heat is released to the container, to the internal air and to the objects immersed in the water (stirrer, electric resistances, thermometer). Since these items have different specific masses and heats, it is better to determine the water equivalent of the calorimeter which is the quantity of water capable of exchanging the same quantity of heat absorbed by the container and the immersed parts. In determining the different calorimetric quantities, this quantity has to be added to the water mass in the calorimeter.

#### *Go ahead as follows:*

- 1. Pour a quantity of water  $m_1$ , for example 250 gr, in the calorimeter.
- 2. Take note of the water temperature  $t_1$  in the calorimeter (20°C, for example).
- 3. Prepare a quantity of water  $m_2$  at a temperature  $t_2$  higher than  $t_1$ .  
E.g.  $m_2 = 250$  g,  $t_2 = 52^\circ\text{C}$ .
- 4. Pour the warm water quickly in the calorimeter and facilitate the heat exchange thanks to the stirrer.
- 5. Take note of the equilibrium temperature  $t$ . E.g.  $t = 34^\circ\text{C}$ .

The heat released by the water mass at a temperature  $t_2$  is  $Q_2 = cm_2(t_2 - t)$   $c$  is the specific heat of water.

The heat absorbed by the water contained in the calorimeter and by the equivalent mass is.

$$Q_1 = c(m_1 + m_e)(t - t_1)$$

Since the heat released is the same of the heat absorbed, equalling the two previous formulas you have the equivalent mass:

$$m_e = \frac{m_2(t_2 - t)}{(t - t_1)} - m_1$$

Considering the data given as an example it will be  $m_e = 71,2$  g.

- 1. Do the measure more than once changing the values of the masses and of the initial temperatures and get the average of the values.

**2) DETERMINING THE SPECIFIC HEAT OF A SOLID*****Go ahead as follows:***

1. Pour a quantity of water  $m_1$  in the calorimeter and measure its temperature  $t_1$ .
2. Measure the mass  $M$  of a solid sample with unknown specific heat and immerse the sample in the water. Porta l'acqua all'ebollizione, controllando la sua temperatura con un termometro.
3. Once the water has reached its boiling point which is indicated as initial temperature of the sample  $t_2$ , quickly immerse the sample in the calorimeter, close and stir.
4. The temperature of the system Water + Sample starts to increase until it stabilizes: take note of the thermal equilibrium value's temperature  $t_e$ .
5. The heat released by the solid to the system is  $Q_C = Mc_x(t_2 - t_e)$ .
6. The heat absorbed by the system, considering the water equivalent, is  $Q_A = c(m + m_e)(t_e - t_1)$ ,  $c$  is the specific heat of water (1 cal/Kgx°C°).
7. Equalling the heat released to the one absorbed, you can determine the specific heat of the solid:

$$c_x = \frac{c(m + m_e)(t_e - t_1)}{M(t_2 - t_e)}$$

**NOTE:** in order to determine the specific heat of a liquid, act following the same steps using an other unknown liquid instead of water and considering the known specific heat of the solid.

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**Optika S.r.l. - Copyright**