Guida didattica / Didactic guide
Guía didáctica

Ver. 1.0.1

Cod. 4209

OPTIKA S.r.l.

Via Rigla, 30 – 24010 Ponteranica (Bergamo) – ITALY
Tel. +39 035 571392 - Fax +39 035 571435

www.optikascience.com

info@optikascience.com

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SPECTROMETER

FEATURES:

- BASE:** made of varnished cast iron. Equipped with a goniometer of 17.5 cm of diameter; divided into 360° with 1° precision. A Vernier calipers allows to measure with an accuracy of 1/10°.
- TELESCOPE:** The telescope is fitted with an achromatic objective of a 178mm focal distance and 15X ocular. With fine focusing. Use the thumbscrew located under the item in order to fix the instrument. A small thumbscrew located on the side of the instrument permits the fine rotation of the telescope along the stage within a small angle.
- COLLIMATOR:** 178mm achromatic objective. Continuously adjustable slit from 0 to 6mm by a side screw.
- PRISM HOLDER STAGE:** continuously adjustable, horizontally and vertically. A side thumbscrew permits the locking of the rotation.
- SUPPLIED ACCESSORIES:** 1 glass equilateral prism; 1 diffraction grid; 1 magnifying lens.

MEASUREMENT OF A PRISM'S REFRACTION INDEX

Firstly put off the grid holder unblocking the screws.

1. Place a monochrome source in front of the collimator's opening. (eg. a mercury-vapor lamp)
2. Observe the image of the opening through the telescope placing the telescope in axis with the collimator, **without the prism**. Focus the image with the ocular of the telescope, evaluate the thickness of the opening's image and if necessary make it thin acting on the side screw.
3. Unblock the telescope and the stage. Perfectly collimate the image of the opening with its cross on the ocular using the rotary motions of the telescope and its fine adjustment. Block the telescope.
4. Rotate slowly the stage in order that the point "zero" of the protractor coincide with the arrow located on the vernier scale.
5. Block the stage.
6. Place the prism as illustrated in figure 1. The opaque face is the one indicated in grey. Respecting this disposition is necessary in order to process the experiment.
7. Unblock the telescope and rotate it slowly and anticlockwise until you can see the image of the opening in the middle of the visual field. Block the telescope and the collimator with the cross using the fine rotation.
8. Unblock the stage and rotate it **slowly** clockwise in a continuous way: The image of the opening moves slowly on the left of the visual field until it reaches *an inversion point*.
9. Note down, from the protractor, the value of the angle reached at the inversion point. This value represents the angle of minimum deviation δ_m .

4209

10. Apply the formula that gives the refraction index:

$$n = \frac{\sin \frac{\Delta + \delta_m}{2}}{\sin \frac{\Delta}{2}}$$

in which with Δ we have indicated the refraction angle of the prism. In this case $\Delta = 60^\circ$ as the prism is equilateral and so

$$\sin \frac{\Delta}{2} = 0,5$$

Example:

Flint prism (nominal refraction index: 1,62)

After having performed all the indicated operations, the result is

$$\delta_m = 48^\circ$$

$$n = \frac{\sin \frac{60 + 48}{2}}{0,5} = 1.618$$

in collusion with the nominal value.


Fig. 1: Refraction index of a prism


Fig. 2: wavelength with the grid

MEASUREMENT OF A WAVELENGTH WITH THE DIFFRACTION GRID

Firstly it is necessary to mount the grid holder support on the stage.

10. Put the diffraction grid in the support.
11. Place the monochrome source in front of the collimator (eg. a mercury-vapor lamp).
12. Place in axis source – collimator – telescope.
13. Block the stage and perfectly collimate the image of the source in the middle of the visual field using the cross and the fine rotation.
14. The position of the image in this disposition represents the maximum of zero order. Note down the value of angle α_0 in this configuration.
15. Keeping the stage blocked, place the diffraction grid so that it is perpendicular to the optical axis and rotate slowly the telescope anticlockwise (or clockwise) until you identify the spectral lines of higher order (figure 2).
16. Note down the corresponding angle α_1 .
17. Calculate the effective angle of the spectral line of higher order, that is $\alpha = \alpha_1 - \alpha_0$.
18. Apply the formula of the grid to calculate the wavelength: $\lambda = p \sin \alpha$, for the maximum of the visible spectral lines.

Practical example: we consider the green line:

Grid with 500 openings / mm; step $p = \frac{1}{500} \text{ mm}$

After having performed the above operations, the result is:

$$\alpha_0 = 324.4^\circ \quad \alpha_1 = 340.5^\circ \quad \alpha = 340.5 - 324.4 = 16.1^\circ$$

$$\lambda = \frac{\sin(16.1)^\circ}{500} \text{ mm} = 554.6 \text{ nm}$$

in collusion with the theoretical value (546.07).

NOTICE

Any differences between the material provided and the drawings showing them is due to technological update

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Optika S.r.l. - Copyright