
5/96-d.R./Sf-

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## Elektroskop Electroscope

Fig. 1

The electrostatic induction experiment is used for basic electrostatic experiments and for indicating DC and AC voltages up to approx. 8 kV.

## 1 Description, technical data, scope of supply

- ① Iron ring for defined field distribution, 18 cm dia.
- ② Support for indicator, insulation-fitted to the iron ring ①
- ③ Indicator, with point bearing, asymmetrically suspended so that its own weight gives the necessary restoring moment
- ④ Insulator with 4-mm socket which is conductively connected with the indicator support ②
- ⑤ Capacitor plate, 4 cm dia., with 4-mm plug pin
- ⑥ Stand pin with 4-mm socket, conductively connected with ring ①

Saddle base (300 11), not included in scope of supply

## 2 Use

Recommended for charging:

Either

High-voltage power

supply unit, 10 kV 521 70

High voltage cable 501 05

Induction plate 542 51

or

Friction rods with rubbing material, e.g.

Friction rods	Rubbing material	Charge polarity
Pair of friction rods PVC and acrylic glass (541 00)	Plastic foil	+; -
Glass tube (541 01)	Newspaper, leather (541 21)	+
Glass tube with bearing pan (541 02)	Newspaper, leather (541 21)	+
Plastic rod (541 04)	Wool or textile fibers	-
Plastic rod with bearing pan (541 05)	Wool or textile fibers	-


Fig. 2

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
Fig. 3

Make sure that insulator ④ is clean and dry. If necessary, clean it with cleaning spirit or alcohol and dry it with a hot air blower.

In case of high humidity or if the apparatus is transferred from a cold into a warm room, dry the electrostatic induction plate and the surroundings with hot air.

Recommended for demonstration experiments: shadow projection using a point-shaped light source (e. g. lamp housing with lamp, 450 60/51)

*Charging an electrostatic induction plate by touching it with a charged body:*

Either

touch the capacitor plate ⑤ on the electrostatic induction plate with rubbed glass or plastic rod,

or (see Fig. 2)

Attach capacitor plate ⑤ to socket a of the high-voltage power supply 10 kV (521 70) and transfer the charge from the capacitor plate to the electrostatic induction plate using the induction plate (542 51).

*Charging by electrostatic induction (see Fig. 3):*

Positive electrostatic induction plate charge: approach the capacitor plate ⑤ of the electrostatic induction plate with rubbed plastic rod, but do not touch the plate. At the same time touch the capacitor plate briefly with your finger. Then withdraw the plastic rod.

For negative electrostatic induction plate charge, use rubbed glass rod.