
7/1989

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The translucent screen can be used as a projection surface for optical phenomena, e. g. for real images, spectra or interference and diffraction patterns. Working with images of weak illumination or in only moderately darkened rooms is greatly facilitated by the use of the screen.

1 Description

Plate of acrylic glass (30 cm x 30 cm x 0.3 cm), one side frosted; on rod (9 cm x 1 cm diam.)

2 Use

The matted surface of the screen is sensitive to grease and perspiration. In order to prevent it being stained, it should, therefore, only be held by the stand rod. By means of the stand rod, it can be fitted in a Leybold boss (301 01) used as a rider for the small optical bench (460 42 or 460 43), into a rider (460 35/45/48/52/55/58) for the optical bench with standard profile (460 32) or the demonstration table according to Prof. Pohl (460 65) or into a stand base (300 01 or 300 02) or saddle base (300 11).

It will depend upon the experimental conditions which side of the translucent screen should be turned towards the source of light. In a completely darkened room (no scattered light, no interfering reflections) the matted side should be turned towards the source of light, and the image should be viewed from the other side in transmitted light.

Whenever a mark shall be made on the translucent screen, it is recommended to cover the screen with a sheet of white paper.

3 Cleaning

Important: Do not use aggressive solvents!

The plate may be cleaned with e. g. "Plexiklar®" made by Messrs. Röhm & Haas, or with a surface-active substance like Tide or Dreft. The surface should be treated with one of these substances using a soft cloth and only moderate pressure. Even if obstinate stains have to be removed, neither harder objects nor specifically fatdissolving agents like benzine or alcohol must be used.