
9/95-Sf-

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Literatur: Zentralkartei Optik (598 641)

The optical disks, designed for vertical assembly in a glancing light beam, permit basic qualitative and quantitative experiments on reflection, refraction and dispersion of light, in conjunction with the clamping sectional models of mirrors, lenses and prismatic and circular bodies.

Literature: Central Card-file IV: "Optics" (598 641, in German)

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#### 1 Safety notes

- Use only water-soluble felt-tip pens for marking and writing on the optical disk.
- Never clean the optical disk and acrylic glass objects with aggressive solvents.  
Recommended: a lint-free cloth moistened with water and a little washing-up liquid.

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#### 2 Description, scope of supply (see Fig. 1)

- ① Optical disk (40 cm x 25 cm), plastic, with imprinted 360°-scale in 1°-divisions; with five 4-mm holes (1.1) for spring clamps ③ - threaded hole (1.2) on rear of panel for stand rod ②
- ② Stand rod (12 cm x 1 cm dia.) with thread for hole (1.2), for mounting the optical disk in stand material
- ③ 2 spring clips with 4-mm pin for sockets (1.1), for mounting the elements ④ through ⑩ on the optical disk
- ④ Mirror holder with planar mirror of glass and spring-mounted metal strip, which can be used both as a concave and a convex mirror
- ⑤ to ⑩ Acrylic glass sectional models, matt on one side
  - ⑤ Right-angled prism
  - ⑥ Trapezoid / planoparallel plate
  - ⑦ Biconvex lens
  - ⑧ Biconcave lens
  - ⑨ Semicircular model
  - ⑩ Circular model
- ⑪ Preformed storage tray (35 cm x 17 cm x 3 cm) for parts ② through ⑩; with cover

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### 3 Operation

Recommended for creating light beam, and for setup:

Lamp 6 V, 30 W (450 51) in lamp housing (450 60) with transformer 6 V, 30 W (e.g. 562 73)

or

Halogen lamp 12 V, 50 W (450 68) in halogen lamp housing 12 V, 50/100 W (450 64) with transformer 12 V, 100 W (e.g. 562 75)

Lens in holder,  $f = 150$  mm (460 08) as condenser

Diaphragm with 5 slits (463 51) for screening out individual colored "light beams"

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Experiment setup with optical disk

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Set up the apparatus as shown in Fig. 2 according to your particular experiment conditions.

When illuminating with the 6 V lamp (Fig. 2), turn insert ① in the lamp housing so that the lamp filament is aligned horizontally (precondition for glancing light along the optical disk);

Set the distance between the lamp and the lens so that you obtain a divergent, parallel or convergent light beam, according to the experiment requirements.

Attach the acrylic glass bodies so that the matt side is facing the optical disk.

To adjust the angles of incidence  $\neq 90^\circ$  (see detail sketches) loosen the screw ② slightly and set the optical disk in the desired position.