
Physics	Chemistry · Biology	Technology



Lehr– und Didaktiksysteme LD Didactic GmbH Leyboldstrasse 1 · D-50354 Huerth

06/05-W97-Sel



1 Description

The pair of resonance tuning forks serves especially for experiments on acoustic beats. It consists of two equal tuning forks on resonance boxes. The frequency of one fork can be changed with a tuning weight.

2 Scope of supply

- 2 Tuning forks on resonance boxes
- 1 Tuning weight
- 1 Striking hammer

3 Technical data


Characteristic frequency:	440 Hz
Tolerance:	approx. 0,5 Hz
Frequency change :	0 Hz (tuning weight: bottom) approx3 Hz (middle) approx20 Hz (top)
Length of tuning fork:	20 cm
Dimensions of resonance boxes:	18 cm \times 9 cm \times 5 cm

Instruction sheet 414 72

Pair of resonance tuning forks (414 72)

- 1 Resonance box
- 2 Tuning fork
- 3 Tuning weight
- 4 Striking hammer

4 Experiments on beats



- Direct the openings of the resonance boxes towards each other.
- Install tuning weight at one tuning fork.
- Strike the two tuning forks one immediately after the other with the striking hammer.

LD Didactic GmbH · Leyboldstrasse 1 · D-50354 Huerth / Germany · Phone (02233) 604-0 · Fax (02233) 604-222 · e-mail: info@ld-didactic.de