Physik

Chemie · Biologie

Technik

LEYBOLD DIDACTIC GMBH

1/1991

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The iris diaphragm is used in optics, above all as a continuously variable stop, to limit a beam of light and screen off scattered light.

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1 Description

The iris diaphragm is inserted into a plastic frame (diam. 120 mm). The diaphragm consists of a circular arrangement of thin metal segments which form an approximately circular hole. By means of a pin on the diaphragm frame the stop aperture can be varied continuously between approx. 1.5 mm and 30 mm of diam.

The stand rod of approx. 10 mm diameter and 80 mm length attached to the frame can be used to clamp the diaphragm to an optical bench, or to stands.

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2 Use


The iris diaphragm should be clamped into the rider of the optical bench with its frame adjusted at right angles to the optic axis. If Leybold bosses (301 01) are used as riders for the small optical bench (460 42 or 43), it should be checked that the stand rod on the iris diaphragm is arranged at the same height as that of the other pieces of equipment. In this manner a coarse adjustment of the iris diaphragm in the optic axis is obtained. Thereafter the width of the opening of the diaphragm should be adjusted to the desired size.

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