8/96-Sf-

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### Small Optical Bench with Short Lateral Bracket Small Optical Bench with Long Lateral Bracket Swivel Joint with Angle Scale

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The optical benches (460 42/43) are used together with Leybold multiclips (301 01) as riders for experiment assemblies for optical experiments requiring adjustment along the axes.

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As well as this, they can always be used as stands wherever an experiment apparatus is required to be particularly mechanically stable, or when devices which can be held in a Leybold multiclip require a guide or need to be moved a defined distance.

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The swivel joint with angle scale (460 40) provides the connection for two small optical benches which must be assembled at a specific and measurable angle to one another for the purpose of carrying out experiments e.g. spectral experiments, experiments on thermal radiation, and for examining diffraction and interference of ultrasonic waves.

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- ① Rectangular-section rail (75 cm long), with stand rod ②, connected so it rotates (360°)
  - (1.1) Scale 0 to 75 cm, with mm-graduation
  - (1.2) Guide groove for screw (4.1) for Leybold multiclips ④ used as riders
- ② Stand rod (10 mm Ø), rotates through 360° with respect to rail  
length of stand rod: 50 cm (for 460 42)  
25 cm (for 460 43)
- ③ Connection between rail ① and stand rod ②, rotatable (360°), fixed using star screw (3.1)
- ④ Leybold multiclips (301 01) used as riders (not included in the scope of supply)


Fig. 1

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### 1.2 Swivel joint with angle scale, see Fig. 2

- ① Clamps, can be rotated by about 190° with respect to one another, with star screws for attaching to two small optical benches  
Clamp (1.1), on the same axis as the indicator ③ and clamp-ing device ④,  
clamp (1.2), on the same axis as the scale ②
- ② Angle scale, -10°...0...+190°, with 1° graduations,  
the scale disk (12 cm dia.) is rotated by rotating the clamp (1.2), or the clamp (1.1); can also be rotated around the axis
- ③ Indicator to show the angle between the optical benches.  
  
The indicator rotates with the clamping device ④ when the clamps (1.1) and (1.2) are rotated; can also be rotated around the axis
- ④ Annular clamping device (dia. 10 mm) for experimental ap-paratus on stand rods with a 10 mm dia.,  
with tensioning screw (4.1) for fixing the apparatus in the clamp,  
rotation of the clamping device together with the indicator ③ when rotating the clamps (1.1) and (1.2)

## 2 Use

Assemble the small optical bench either in a large stand base (300 01), or in a bench clamp (301 06 or 07), (see Fig. 3.1-3.3).

Always mount the Leybold-multiclamps (301 01) or clamps ① of the swivel joint (460 40) as shown in Fig. 4.1/4.2 so that the screws engage in the guide groove (1.2) of the optical bench.

Mount the optical elements with 10 mm dia. stand rod, e.g. lamp housing (450 60), 12 volt halogen lamp, 50 W, 100 W (450 64), lenses (460 01-12), adjustable joint (460 14), iris diaphragm (460 24), etc. - to make approximate adjustments to the axis, clamp the arrangement in such a way that end surface of the stand rod is flush with the underside of the clamp (see Fig. 4.1).

Before connecting the optical benches using the swivel joint with the angle scale (460 40), align them horizontally and bring them to the same height - see Fig. 5.1/2 .


Fig. 3.1


Fig. 3.2


Fig. 3.3

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**Correct mounting**

- of a Leybold multiclamp (303 01) for supporting optical elements on rods (Fig. 4.1)
- of a swivel joint with angle scale (460 40) for connecting two optical benches (Fig. 4.2)


Fig. 5.1

Fig. 5.1/5.2

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Connection of two small optical benches (460 42 or 43) using a swivel joint with angle scale (460 40)  
 Connection of the swivel joint to two optical benches mounted at the same height and aligned horizontally as shown in Fig. 5.1


Fig. 5.2