
Gebrauchsanweisung Directions for Use Mode d'Emploi

46312

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The convexo-concave mirror is suitable for use in a number of optical experiments, and is especially useful for demonstrating the path of rays and the formation of the image by means of spherical mirrors. Again, it is applicable in those optical arrangements in which the infrared rays absorbed are to be kept to the minimum.

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1. Description

Both mirrors are mounted in a frame in such a way that the convex mirror is visible from one side and the concave mirror from the other. The mirrors are approx. 70 mm in dia. with focal lengths of approx. ± 250 mm. The concave mirror is front-silvered.

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The stand rod (approx. 10 mm ϕ , 80 mm long), which is fitted to the frame of the lens, is used for fixing on to an optical bench or on to stands.

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2. Use

Clamp mirror the frame of which is square with the optical axis in the rider of the optical bench. In case Leybold bosses (301 01) are used as riders for the small optical bench (460 42 or 460 43), it is advisable to clamp the stand rod of the mirror as well as those of the other apparatus equally high. Fit all bosses on the same side of the optical bench. This provides a rough adjustment of the mirror in the optical axis.

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