*Guida didattica / Teaching guide*  
*Guía didáctica*

Ver. 1.0.1

---

**Cod. 3006**

---

**OPTIKA S.R.L.**

VIA RIGLA, 30 – 24010 PONTERANICA (BERGAMO) – ITALY  
Tel. +39 035 571392 - Fax +39 035 571435

[www.optikascience.com](http://www.optikascience.com)

[info@optikascience.com](mailto:info@optikascience.com)

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
## DEVICE FOR STUDYING WAVES

Waves are an essential part of our daily life. sound and light waves are deciphered by our ears and eyes: they are the pressure waves and the electromagnetic waves, respectively. we normally meet mechanic vibration waves, radio waves, x-rays, infrared and ultraviolet, the guided waves propagation and the delay lines. at a basic level, all particles properties can be described by schrodinger's probability waves. The device for transverse waves, provides a example of transverse waves' propagation, where the angle of each vibrant rod has properties related to the wave.

### *THE MACHINE FOR TRANSVERSAL WAVES*


Picture 1

The wooden machine for transversal waves is made by inserting a nylon thread in one of the two holes located in each wooden cylindrical rod and then through one of the plastic rods (2.5 cm length). This must be repeated for all 80 wooden cylindrical rods.

Let loose about 20 cm of the thread in order to make a knot, then put the thread into the second hole of all wooden rods, which have plastic spacers before the second hole. Then cut the thread and tie its extremities so as to make another knot. This must be done by joining the thread extremities and making at least three coincident knots so as to prevent the rods from moving when the thread is tensioned.

The handles are fixed by means of these knots as shown in picture 2.

Picture 2

It can be prepared also by fixing an extremity to a wall while the other is kept in one hand, keeping it tensioned and in horizontal position. Otherwise, two persons may keep the thread tensioned and horizontal by pulling both handles. Once it has been assembled and stabilized, all rods should point to the same direction that can be vertical, horizontal or any other.

Each rod's inertia momentum can be increased by tightening the nuts at the wooden rods' extremities.

It may be necessary to wrap up the wooden rod by means of an adhesive tape in order to support the nut.

### **NOTICE**

Small differences in the provided components and their drawings are justified by the technological updates.

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Optika S.r.l. - Copyright